

Fitting & Care Instructions

Model: DCT-4351



Instructions on putting on brace:



Step 1: Center brace on side of knee and feed the bottom calf strap through the D-Ring and tighten. If strap is too short there is an extension piece that is provided with the brace



Step 2: Feed upper thigh strap through the D-Ring and tighten. If strap is too short there is an extension piece that is provided with the brace



Step 3: Feed the bottom strap around the back of the calf and pull through D-Ring. If the strap is too long please revert to opposite side to shorten strap



Step 4: Feed the top strap around the thigh and pull through D-Ring. If the strap is too long please revert to opposite side to shorten



Step 5: Once all straps are tightened, use provided wrench to turn dial counter clockwise. Do small 1/4 inch turns and walk to check relief. Repeat until the preferred amount of unloading is reached. Adjust top and bottom to desired comfort.



Step 6: Once all straps are tightened, use provided wrench to turn dial counter clockwise. Do small 1/4 inch turns and walk to check relief. Repeat until the preferred amount of unloading is reached. Adjust top and bottom to desired comfort.

See reverse side for strap fitting adjustment instructions.

Strap Adjustment Instructions:



Step 1: After initial fitting if any of the straps were too long, lay the brace out open on a table



Step 2: Locate the straps that were initially too long. On the inside you can adjust the length of the strap

Adjust the strap length by pulling the strap toward the middle.

Step 3: If the straps are too short attach the provided extension panel to the end of each strap.



CARE INSTRUCTIONS

Hand wash with mild detergent and air dry.

WARRANTY

The DCT-4351 is covered by a 6 month warranty against all mechanical or manufacturer defects, from the date of purchase.