

Trend Flex Wrist & Thumb

DCT-WTFLEX The Trend Wrist & Thumb Flex is designed for post-operative rehabilitation, tendon and ligament injuries, thumb stiffness, and contracture management, this brace allows adjustable thumb flexion and extension for controlled or progressive positioning. It can also be used as a traditional thumb spica for stabilization and immobilization.



Customizable Thumb Positioning

An integrated adjustment system allows the thumb position to be precisely chosen through controlled flexion and extension. With the provided adjustment tool, the thumb component can be positioned to the desired setting for individualized support.



Highlights

- Clinician-adjustable thumb positioning
- Supports controlled flexion, extension, and/or abduction
- Adjustable wrist and thumb stabilization
- Locking mechanism for customized positioning
- Supports progressive positioning

Indications

- Thumb CMC osteoarthritis
- Wrist and thumb sprains or strains
- Thumb CMC/MCP instability and ligament injuries
- Post-operative or post-traumatic thumb rehabilitation
- De Quervain's tenosynovitis

PDAC Approved: L3915/L3916